

SmartPractice[®].com **Glove Sizing Guide**

Instructions:

Size is indicated by the width of your hand just below your first knuckles.


Place your right hand flat (palm down) on the guide with your fingers together and thumb spread. Your index finger should be flush against the heavy line, but not over it.

PRINT OUT CHART TO USE

(cannot be used when on your screen)

When printing out this Glove Sizing Chart, please be sure to **deselect** either the “fit to page” or “shrink to fit” option within the print menu.

*Sizing for SmartPractice[®]
Brand Gloves*



To ensure accuracy: make sure that this width measures $3 \frac{15}{16}$ " or 10 cm as shown.